

# CUTERA®

A New Energy in Aesthetics

## PORTFOLIO OF PRODUCTS

Cutera is a leading provider of laser and other energy-based aesthetic systems for practitioners worldwide. Since 1998, Cutera has been developing innovative, easy-to-use products that enable physicians and other qualified practitioners to offer safe and effective aesthetic treatments to their patients.




# CUTERA PORTFOLIO OF PRODUCTS



## secret™ PRO

### TAKING ANTI-AGING TO NEW DEPTHS

Secret PRO provides skin revitalisation by uniquely combining two clinically proven technologies in a single device. Using fractional CO<sub>2</sub> for skin resurfacing and radiofrequency microneedling for deep dermal remodeling, Secret PRO effectively treats from the epidermis to the deeper dermal layers independently or together.<sup>6-8</sup>



SCAN TO LEARN MORE



Photos courtesy of Somenek + Pittman MD. Individual results vary. Images are for illustrative purposes only. A consultation is required to determine suitability. vary.

For device usage and safety guide, visit [cutera.com.au](http://cutera.com.au)

## secret™ DUO

### DYNAMIC FACIAL RESURFACING AND REMODELING

Secret DUO pairs RF microneedling with a non-ablative 1540 nm laser in a single platform to target the common signs of facial aging such as fine lines, texture, striae, and scars. Each modality, used individually or in combination, offers customisable treatments and optimal flexibility to address a variety of skin types and conditions.<sup>6-9</sup>



SCAN TO LEARN MORE



Photos Courtesy of Evren Gökeşme, MD. Individual results vary. Images are for illustrative purposes only. A consultation is required to determine suitability.



For device usage and safety guide, visit [cutera.com.au](http://cutera.com.au)

## secret™ RF

### SKIN REMODELING FROM THE INSIDE OUT

Secret RF is a tailored approach to skin revitalisation combining bipolar RF energy and microneedling with both semi-insulated and non-insulated needles to treat all skin types. Address the common signs of aging via hemostasis and coagulation as well as large pores, scars, stretch marks, and photoaging.<sup>6-8</sup>



SCAN TO LEARN MORE



Photos courtesy of J. Hsu, MD. Individual results vary. Images are for illustrative purposes only. A consultation is required to determine suitability.

For device usage and safety guide, visit [cutera.com.au](http://cutera.com.au)

# CUTERA PORTFOLIO OF PRODUCTS

## xeo+

### UNLOCK LIMITLESS POSSIBILITIES

Transform your aesthetic practice with xeo+, the multi-platform solution designed to deliver more. From skin rejuvenation and hair removal to vascular and pigmented lesion treatments, xeo+ offers unmatched versatility and cutting-edge precision.



#### COMPREHENSIVE TREATMENTS

Address diverse patient needs with one powerful platform.



#### INCLUSIVE CARE

Deliver trusted treatments, safe for all skin types.



#### NON-INVASIVE PRECISION

Achieve remarkable results without downtime.



#### ADAPTABLE GROWTH

Expand and customise as your practice advances.



SCAN TO  
LEARN MORE

Features Cutera's  
Signature  
LASER Genesis®  
Treatment



Model, not actual patient.

For device usage and safety guide, visit [cutera.com.au](http://cutera.com.au)



Photos courtesy of Nina Desai, DO. Individual results vary. Images are for illustrative purposes only. A consultation is required to determine suitability.



Photos courtesy of Terri Raymond, PA-C. Individual results vary. Images are for illustrative purposes only. A consultation is required to determine suitability.

## enLIGHTen®SR

### SKIN REVITALISATION RE-IMAGINED

A novel picosecond platform featuring two proven technologies, 532 nm and 1064 nm, and a high pulse energy for safe and effective treatments for unwanted pigmentation and skin revitalisation concerns. Expand applications to include scarring and difficult-to-treat skin concerns with PICO Genesis FX.



Photos courtesy of S. Ronan, M.D. Individual results vary. Images are for illustrative purposes only. A consultation is required to determine suitability.

SCAN TO  
LEARN MORE



Model, not actual patient.

For device usage and safety guide, visit [cutera.com.au](http://cutera.com.au)

# CUTERA PORTFOLIO OF PRODUCTS



## THE FIRST AND ORIGINAL TGA-CLEARED ENERGY DEVICE FOR THE TREATMENT OF MILD TO SEVERE ACNE

AviClear is the first and original TGA-cleared energy device for the long-term treatment of mild, moderate, and severe acne.\* Using a 1726 nm wavelength, AviClear targets acne at its source, and can significantly reduce acne in just three 30-minute treatment sessions.<sup>1-3</sup>



SCAN TO LEARN MORE



Photos courtesy of Cutera Inc. Individual results vary. Images are for illustrative purposes only. A consultation is required to determine suitability.

For device usage and safety guide, visit [cutera.com.au](http://cutera.com.au)

## excel<sup>®</sup>V+

### PRECISION & PERFORMANCE IN VASCULAR AND PIGMENTATION TREATMENTS

excel V+ is the latest-generation laser technology for vascular and pigmentation treatments. With design input from leading dermatologists, excel V+ delivers the power, precision, and performance to safely and effectively treat indications, from challenging vascular and pigmentary conditions to today's most common skin concerns.<sup>10</sup>



SCAN TO LEARN MORE



Photos courtesy of E.V. Ross, M.D. Individual results vary. Images are for illustrative purposes only. A consultation is required to determine suitability.



For device usage and safety guide, visit [cutera.com.au](http://cutera.com.au)



## enLIGHTen<sup>®</sup>

### POWERFULLY VERSATILE PICO TECHNOLOGY

A versatile portfolio of picosecond technology featuring multiple wavelengths, including the first true-red picosecond wavelength, and industry-leading specifications to address multicolour tattoos, the treatment of benign pigmented lesions and PICO Genesis, and PICO Genesis FX skin revitalisation.<sup>13-15</sup>



SCAN TO LEARN MORE



Photos courtesy of E.V. Ross, M.D. Individual results vary. Images are for illustrative purposes only. A consultation is required to determine suitability.

For device usage and safety guide, visit [cutera.com.au](http://cutera.com.au)

# CUTERA PORTFOLIO OF PRODUCTS



## 'truFlex™

### REDEFINING MUSCLE SCULPTING

Tone and firm muscle with truFlex, a customisable muscle stimulation device featuring multidirectional stimulation (MDS) technology to strengthen, tone, and define muscles on the abdomen, obliques, glutes, and thighs. Treat up to eight areas at once and increase muscle mass by up to 30%.<sup>4</sup>



SCAN TO LEARN MORE



Photos courtesy of M. Kaga, M.D. Individual results vary. Images are for illustrative purposes only. A consultation is required to determine suitability.

For device usage and safety guide, visit [cutera.com.au](http://cutera.com.au)

## 'truSculpt®

### SCULPT YOUR BEST BODY

Target stubborn fat and give your patients personalised results with one of the fastest body sculpting treatments on the market. With an average 24% fat reduction per treatment area, this award-winning technology utilises heat through monopolar radiofrequency to break down fat so you can deliver results with high patient satisfaction.<sup>5</sup>



SCAN TO LEARN MORE



Photos courtesy of N. Desai. Individual results vary. Images are for illustrative purposes only. A consultation is required to determine suitability.



For device usage and safety guide, visit [cutera.com.au](http://cutera.com.au)



## 'truBody®

### REMOVE FAT, REBUILD MUSCLE

truBody is the only customisable fat-reducing and muscle-toning aesthetic duo that can accomplish results with no downtime in 30 minutes or less. truBody combines truFlex and truSculpt to deliver a customised, convenient, and complete solution that addresses today's in-demand treatments and delivers results with high patient satisfaction.



SCAN TO LEARN MORE



Photos courtesy of M. Kaga, M.D. Individual results vary. Images are for illustrative purposes only. A consultation is required to determine suitability.

For device usage and safety guide, visit [cutera.com.au](http://cutera.com.au)

# About Cutera

Cutera is a leading provider of aesthetic and dermatology solutions for practitioners worldwide. For over 25 years, Cutera has strived to improve lives through medical aesthetic technologies that are driven by science and powered through partnerships.

For more information, call 1800 CUTERA (288 372) or visit [cutera.com.au](http://cutera.com.au)



Over 25 years, we've established a legacy of ingenuity, engineering, and proven partnership that's brought new energy to practices around the world.

## CUTERA®

Cutera ANZ

24/1 O'Connell Street, Sydney NSW 2000, Australia

1800 CUTERA (288 372) | [aumarketing@cutera.com](mailto:aumarketing@cutera.com) | [cutera.com.au](http://cutera.com.au)

### References:

1. Data on file. FDA clearance study. Cutera Inc.
2. O'Neill AM, Gallo RL. Host-microbiome interactions and recent progress into understanding the biology of acne vulgaris. *Microbiome*. 2018;6:177.
3. Sakamoto FH, et al. Selective photothermolysis to target sebaceous glands: theoretical estimation of parameters and preliminary results using a free electron laser. *Lasers Surg Med*. 2012;44(2):175-183.
4. Ronan SJ, A novel bio-electric current stimulation device for improvement of muscle tone. 7/2019.
5. Taub, A, Bartholomeusz J., Ultrasound Evaluation of a Single Treatment with a Temperature Controlled Monopolar Radio Frequency Device for the Improvement of Localised Adiposity on the Abdomen and Flanks. *J Drugs Dermatol*. 2020;19(1):28-34.
6. Tierney E and Hanke W. Fractionated Carbon Dioxide Laser Treatment of Photoaging: Prospective Study in 45 Patients and Review of the Literature. *Dermatol Surg*. 2011;37(9):1279-1290.
7. Yadav S and Singh A. Microneedling: Advances and Widening Horizons. *Indian Dermatol Online J*. 2016;7(4):244.
8. Ryu H-W, et al. Clinical Improvement of Striae Distensae in Korean Patients Using a Combination of Fractionated Microneedle Radiofrequency and Fractional Carbon Dioxide Laser. *Dermatol Surg*. 2013;39(10):1452-1458.
9. Cho SI, et al. Evaluation of the Clinical Efficacy of Fractional Radiofrequency Microneedle Treatment in Acne Scars and Large Facial Pores. *Dermatol Surg*. 2012;38(7):1017-1024.
10. Data on File. Cutera, Inc. Brisbane, CA. 2016. (510(k) clearance, K153671). 2. excel V+ Operator Manual. Cutera, Inc. D2325, Rev C, 12/21.
11. Xeo Operator Manual. Cutera, Inc. D0075, Rev L, 02/12.
12. Data on File. Cutera, Inc. Brisbane, CA. 2016. (510(k) clearance, K153671).
13. Data on File. Cutera, Inc. Brisbane, CA. 2018. (510(k) clearance, K182997).
14. Data on File. Cutera, Inc.
15. enlighten III Treatment Guidelines. Cutera, Inc. D1967, Rev F, 06/18. 4. Bernstein EF. Laser Tattoo Removal. *Semin Plast Surg*. 2007;21(3):175-192.
16. excel HR Operator Manual. Cutera, Inc. D1492, Rev C, 11/14.

Important safety and risk information about Cutera products are available by contacting Cutera at [www.cutera.com/us-en/global-offices](http://www.cutera.com/us-en/global-offices).

©2025 Cutera ANZ. All rights reserved. APO04141 rA (05/25) .